

Materials

RNeasy Mini Kit (Qiagen)
Qias shredder columns (Qiagen)
MMLV (NEB)
Exonuclease I (NEB)
Rnasin plus (Promega)
SYBR Green Real time PCR mastermix without ROX (MCLab)
ddH₂O
200 mM MOPS pH 7.5
500mM EDTA
5M NaCl
Brij-700 (5% solution in water)
SYBR Green dye (Invitrogen)
dT₁₆ primer
V, J and C oligonucleotides (see spreadsheet)

Equipment

Real time PCR machine capable of multiple reads during annealing step

We recommend the MJ Research Opticon [Bio-Rad] but other machines such as the ABI7700 may be used as well.

Consumables for Opticon:

Optical caps
Strip tubes (white)
96 well plate (white)

Optional:

- a. Nanodrop spectrophotometer or Agilent Bioanalyzer
- b. Regular PCR machine for RT and non-real time amplification

Method

1. Start with lymphocyte samples (10^5 to 10^7 cells), obtained via
 - a. Dissociated spleen or thymus
 - b. Ficoll prepped PBMCs from blood or a buffy coat
 - c. FACS or Miltenyi sort
2. Pellet cells in eppendorf tubes by spinning 10 minutes at 3000 rpm in refrigerated microfuge. Aspirate supernatant
3. Disrupt pellet in buffer RLT (Qiagen). Use 350 μ L if less than 5×10^6 cells in sample, 600 μ L if 5×10^6 cells or more in sample. Place RLT suspension through Qias shredder column (1 minute spin at top speed in microfuge) to homogenize. Samples may be frozen at -80°C at this point.
4. Prepare RNA using Rneasy columns following the Qiagen protocol. RNA should be eluted with 34 μ L RNase-free water; final volume will be about 30 μ L. Expected yield will be at least 1mg/ 10^6 cells. Freeze RNA at -80°C .
5. Optional ways to check RNA yields:

- a. Nanodrop spectrophotometer
 - b. Oligogreen assay (Molecular Probes) according to manufacturer's directions. Use a 100 μ l assay volume and read FAM fluorescence with the real time PCR machine.
 - c. Agilent bioanalyzer chip.
6. Combine 14 μ l RNA with 1 μ l dT18 (15uM) and 1 μ l Rnasin plus (40 units/ μ l; Promega) in a 200 μ L strip tube. [Note: regular Rnasin is not heat-resistant and must be added in step 7 instead.] Heat to 65°C for 10 minutes in PCR machine then remove promptly and chill on wet ice.
 7. To each tube add:
 - 2 μ l 10x MMLV buffer
 - 1 μ l 10mM dNTPs
 - 1 μ l MMLV RT (200 units/ μ L; NEB)
 8. Incubate 37°C x 60 minutes, then 95°C x 5 minutes. Freeze leftover RTs at -80°C.
 9. First round V-C amplification
 - 50 μ l mastermix (MCLab)
 - 2 μ l V-specific primer (15 μ M)
 - 2 μ l universal C primer (15 μ M)
 - 4 μ l RT reaction
 - 42 μ l water
 - 100 μ L total reaction
 Program:
 - 95°C 2 minutes
 - 35 cycles of
 - 95°C 15 seconds
 - 60°C 30 seconds
 - final extension
 - 72°C 10 minutes
 Freeze leftover VC reactions at -20°C.
 10. Second round V-J amplification
 - Set up this reaction for each J gene required:
 - 25 μ l mastermix (MCLab)
 - 1 μ l V-specific primer (15uM)
 - 1 μ l J primer (15uM)
 - 0.5 μ l V-C PCR reaction
 - 22.5 μ l water
 - 50 μ L total reaction
 Program:
 - 95°C 5 minutes
 - 10 cycles of
 - 95°C 15 seconds
 - 60°C 30 seconds
 - final extension

72°C 10 minutes

Follow fluorescence to make sure reactions did not reach saturation.

If desired, clonal standards may also be amplified from plasmids, in order to create standards for the annealing reactions.

9. Remove single stranded DNA from reaction (excess primers, asymmetric PCR products, cDNA remnants):

Add 1 μ l Exonuclease I (NEB) to each tube.

Incubate at 37°C for 60 minutes.

10. To each well add 12.5 μ l of reannealing buffer. Mix well; spin plate if there are bubbles. Protect from light as possible.

For 1 mL annealing buffer:

552 μ l 200mM MOPS pH 7.5

110 μ l 500mM EDTA

276 μ l 5M NaCl

34 μ l 5% Brij-700 (protect stock from light)

28 μ l dye concentrate (10,000x SYBR stock from Molecular Probes/Invitrogen, diluted 1:10 in DMSO; protect from light and store at -20°C)

Final concentrations in annealing reaction:

20mM MOPS

10mM EDTA

250mM NaCl

0.03% Brij-700

5x SYBR

11. Check samples by doing 70-99°C melting curve on 10 μ l aliquots
12. Perform annealing reaction, at $T_m - 3.5^\circ\text{C}$ (just below edge of peak on negative first derivative plot of melting curve). Split reaction into 25 μ l sample and 25 μ l reference on either side of a 96 well plate. Up to 32 samples may be measured at once, using the left 4 columns and the right 4 columns of the plate. Fluorescence reads should be every 5 seconds to every 20 seconds, and annealing time can be 1-6 hours depending on expected diversity of sample.
 - a. Make sure that the gradient is set so the rightmost reference well will stay below the annealing temperature.
 - b. Make sure the machine's temperature regulation is set to block and not sample.
 - c. Do melting curves at conclusion of annealing step.